



**KRISHNA VISHWA VIDYAPEETH  
(DEEMED TO BE UNIVERSITY), KARAD**

**UNIVERSITY EDUCATION AND TECHNOLOGY CELL &  
CENTRE OF DISTANCE AND ONLINE EDUCATION KVV**

# **Assay of Ammonium Chloride**

**By**

**Ms. Monali M Shewale**

**Assistant Professor,**

**Krishna Institute of Pharmacy, Karad**

# Objectives


Objectives of this topics is :

1. To determine the percentage purity of given sample of ammonium chloride.

**Titration** :Titration is defined as method of analysis involved in determination of volume of test Solution (unknown concentration solution) by using standard solution (Known conc solution).

**Titrant**: The Solution of Known concentration which is usually taken in a burette.

**Titrate** : The solution of whose concentration is to be determined (unknown Concentration) & usually taken Conical flask



- **Indicators:** A reagent is used to indicate the end point of titration by bring out clear visual colour change of Solution on addition of little amount.



## End point :

The point in the titration at which reaction between two solution is just completed at which the indicator can show sharp colour change is called as end point.



- **Aim:**

To perform assay of given sample of ammonium chloride IP.

- **Requirements:**

- 1. Apparatus:**

Beaker, conical flask, pipette, burette, measuring cylinder, burette stand and holder, glass funnel.

- 2. Chemicals:**

Ammonium chloride, sodium hydroxide, formaldehyde solution, phenolphthalein indicator solution.


- **Principle:**

- 1) The assay of ammonium chloride is carried out by aqueous acid-base titration which based on neutralization reaction.
- 2) It is a direct type of titration.
- 3) Ammonium chloride is treated with formaldehyde solution to liberate equivalent amount of hydrochloric acid.
- 4) The liberated hydrochloric acid is subsequently titrated with 0.1 M sodium hydroxide solution.

- **Theory:**

- (A) Ammonium chloride (Analyte):

- 1. Molecular formula:  $\text{NH}_4\text{Cl}$
- 2. Molecular weight: 53.53.
- 3. Melting point:  $338^\circ\text{C}$
- 4. Description: Ammonium chloride is a white crystalline salt that is highly soluble in water. Solutions of ammonium chloride are mildly acidic.
- 5. Pharmacological activity:
  - (i) Ammonium chloride is used as an expectorant in cough medicine.
  - (ii) Ammonium salts are irritant to the gastric mucosa and may induce nausea and vomiting.





- **Reagent preparation :**
- **(a) 0.1 M Sodium hydroxide:**
- It may be prepared by dissolving 4 g of sodium hydroxide in sufficient carbon dioxide-free water to produce 1000 mL.




- **Assay Procedure:**

- 1. Weigh accurately about 0.1 g of ammonium chloride.
- 2. Dissolve in 20 mL of water.
- 3. Add a mixture of 5 mL of formaldehyde solution, add 20 mL of water.
- 4. Add few drops of phenolphthalein indicator.
- 5. After few minutes, titrate slowly with 0.1 M sodium hydroxide until faint pink colour appears.

- **Reaction:**



Urotropine



Observation table :

<u>Sr.no</u>	Burette Reading (ml) of <u>NaOH</u>		
	Initial	Final	Mean
1.	0.0 ml	12 ml	$\frac{12+15+13}{3}$ $= 13 \text{ ml}$
2.	12 ml	27 ml	
3.	27 ml	38 ml	

# Calculation :

- **Step I: (Amount of ammonium chloride)**

1 mL of 0.1 M NaOH = 0.005349 g of ammonium chloride

13mL(Burette reading of assay) of 0.1 M NaOH = A g of  $\text{NH}_4\text{Cl}$

$$\text{Ag} = \frac{13\text{mL} \times 0.1 \times (0.005349)}{1 \times 0.1}$$
$$= \quad \mathbf{0.06 \text{ gm}}$$

- **Step II: (% Purity)**

For 0.1 g of ammonium chloride 100%

A g of ammonium chloride  $\rightarrow$  p%

$$\text{Hence,} \quad p = \frac{(100) \times 0.06}{0.1}$$
$$= \quad \mathbf{60\%}$$



- **Limit:**
- Ammonium chloride contains not less than 99.0% and not more than 100.5% of  $\text{NH}_4\text{Cl}$ , calculated on the dried basis.



- 1. The percent purity of given sample of ammonium chloride was found to be **60%**.



- In today's class we all have completed with Aim , Principal, Theory , procedure , Observation table and calculations of assay.
- From this practical students can able to find out % purity of any sample.

Thanks for support and guidance:

**Prof. (Dr) Namdeo R. Jadhav**

Dean, Krishna Institute of Pharmacy,  
Krishna Vishwa Vidyapeeth, Karad

**Dr. M. V. Ghorpade**

Registrar,  
Krishna Vishwa Vidyapeeth, Karad

And above all

**Hon. Dr. Sureshji Bhosale**

Chancellor, KVV,

Chairman, Krishna Charitable Trust, Karad



# **END CREDITS**

**Concept, script and cast:**

**Location:**

**Graphics:**

**Narration:**

**Editor:**